

## Product Datasheet

# PKA II $\beta$ reg (phospho Ser113) antibody (orb1415230)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to PKA II $\beta$ reg (phospho Ser113).   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IHC-P, WB  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <b>P31323</b>  |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>Immunohistochemistry: 1/100 - 1/300.<br>ELISA: 1/10000. Not yet tested in other applications.     |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |