

## Product Datasheet

### RUNX1 (phospho Ser249) antibody (orb1415266)

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RUNX1 (phospho Ser249).                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Tested Applications</b> | WB                                                                                                                                     |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Uniprot ID</b>          | <b>Q01196</b>                                                                                                                          |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>ELISA: 1/5000. Not yet tested in other applications.                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                        |