

## Product Datasheet

# PI 3-kinase p85/p55 (phospho Tyr467/199) antibody (orb1415313)

<b>Description</b>	Rabbit polyclonal antibody to PI 3-kinase p85/p55 (phospho Tyr467/199).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC-P, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P27986</a> , <a href="#">Q92569</a>
<b>Dilution Range</b>	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.