

Product Datasheet

p53 (phospho Thr18) antibody (orb1415323)

Description	Rabbit polyclonal antibody to p53 (phospho Thr18).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04637
Dilution Range	IF: IF: 1:50-200 WB 1:500-2000, IHC 1:50-300 IHC 1:50-300
Expiration Date	12 months from date of receipt.