

Product Datasheet

NFκB-p105 (phospho Ser907) antibody (orb1415348)

Description	Rabbit polyclonal antibody to NFκB-p105 (phospho Ser907).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P19838
Dilution Range	WB: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.