

## Product Datasheet

### FoxO1 (phospho Ser256) antibody (orb1415408)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to FoxO1 (phospho Ser256).   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IF, IHC-P, WB   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                    |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q12778</a>  |
| <b>Dilution Range</b>      | WB: Western Blot: 1/500 - 1/2000.<br>Immunohistochemistry: 1/100 - 1/300.<br>Immunofluorescence: 1/200 - 1/1000.<br>ELISA: 1/20000. Not yet tested in other applications. |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |