

Product Datasheet

HSPB8/HSP22 antibody (orb1415672)

Description	Mouse monoclonal antibody to HSPB8/HSP22.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9UJY1
Dilution Range	WB: WB 1:1000-2000
Expiration Date	12 months from date of receipt.