

Product Datasheet

PEPC (Phospho Ser15) antibody (orb1415792)

Description	Rabbit polyclonal antibody to PEPC (Phospho Ser15).
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P20142
Dilution Range	WB: WB: 1:1000-2000
Expiration Date	12 months from date of receipt.