

Product Datasheet

Histone H4 (Phospho Ser1) antibody (orb1415866)

Description	Rabbit polyclonal antibody to Histone H4 (Phospho Ser1).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P62805
Dilution Range	WB: WB: 1:1000-2000
Expiration Date	12 months from date of receipt.