

Product Datasheet

Histone H3 (Phospho Thr32) antibody (orb1415876)

Description	Rabbit polyclonal antibody to Histone H3 (Phospho Thr32).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P84243 , Q71DI3 , P68431
Dilution Range	WB: WB: 1:1000-2000
Expiration Date	12 months from date of receipt.