

## Product Datasheet

# Histone H3 (Phospho Thr11) antibody (orb1415878)

<b>Description</b>	Rabbit polyclonal antibody to Histone H3 (Phospho Thr11).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P84243</a> , <a href="#">Q71DI3</a> , <a href="#">P68431</a>
<b>Dilution Range</b>	WB: WB: 1:500-1000
<b>Expiration Date</b>	12 months from date of receipt.