

Product Datasheet

GATA-2/3 (Acetyl-Lys336/304) antibody (orb1416003)

Description	Rabbit polyclonal antibody to GATA-2/3 (Acetyl-Lys336/304).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P23771 , P23769
Dilution Range	WB: WB: 1:500-10000 ELISA: 1:10000
Expiration Date	12 months from date of receipt.