

Product Datasheet

E2F-3 (Acetyl-Lys168) antibody (orb1416004)

Description	Rabbit polyclonal antibody to E2F-3 (Acetyl-Lys168).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	000716
Dilution Range	WB: WB: 1:500-10000 ELISA: 1:10000
Expiration Date	12 months from date of receipt.