

## Product Datasheet

# Cytokeratin-pan (Acetyl Lys194) antibody (orb1416036)

<b>Description</b>	Rabbit polyclonal antibody to Cytokeratin-pan (Acetyl Lys194).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P13647, Q86Y46, Q7RTS7, P12035, Q3SY84, P08729, P35908, P48668, Q14CN4, P04259, P02538, Q01546, O95678, Q5XKE5</b>
<b>Dilution Range</b>	WB: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>Expiration Date</b>	12 months from date of receipt.