

Product Datasheet

Ku-70 (Acetyl Lys542) antibody (orb1416045)

Description	Rabbit polyclonal antibody to Ku-70 (Acetyl Lys542).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P12956
Dilution Range	WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Expiration Date	12 months from date of receipt.