

Product Datasheet

Ref-1 (Acetyl Lys6) antibody (orb1416056)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Ref-1 (Acetyl Lys6). |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P27695 |
| Dilution Range | WB: Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications. |
| Expiration Date | 12 months from date of receipt. |