

Product Datasheet

Cleaved-Notch 1 (V1754) antibody (orb1416118)

Description	Rabbit polyclonal antibody to Cleaved-Notch 1 (V1754).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P46531
Dilution Range	WB: WB 1:500-2000, IHC-p 1:50-300, IF 1:50-300
Expiration Date	12 months from date of receipt.