

Product Datasheet

Hirudin antibody (orb1416361)

Description	Rabbit polyclonal antibody to Hirudin.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01050
Expiration Date	12 months from date of receipt.