

Product Datasheet

EBV Nuclear Antigen antibody (orb1417079)

Description	Rabbit polyclonal antibody to EBV Nuclear Antigen.
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q1HVF7
Expiration Date	12 months from date of receipt.