

Product Datasheet

MEX3A antibody (orb1418370)

Description	Rabbit polyclonal antibody to MEX3A.
Species/Host	Rabbit
Reactivity	Bovine, Canine, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	A1L020
Expiration Date	12 months from date of receipt.