

Product Datasheet

Otoferlin antibody (orb1418897)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Otoferlin. |
| Species/Host | Rabbit |
| Reactivity | Bovine, Canine, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC-Fr, IHC-P |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | Q9HC10 |
| Expiration Date | 12 months from date of receipt. |