

Product Datasheet

Simian Rotavirus VP4 antibody (orb1418918)

Description	Rabbit polyclonal antibody to Simian Rotavirus VP4.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P12473
Expiration Date	12 months from date of receipt.