

## Product Datasheet

### Influenza A H1N1 Neuraminidase antibody (orb1419294)

|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Influenza A H1N1 Neuraminidase.                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Tested Applications</b> | ICC, IF, IHC-Fr, IHC-P, WB                                                                                                             |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                        |