

Product Datasheet

Hepatitis C Virus 1b Core protein p19 antibody (orb1419356)

Description	Rabbit polyclonal antibody to Hepatitis C Virus 1b Core protein p19.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P26662
Expiration Date	12 months from date of receipt.