

Product Datasheet

GALNTL4 antibody (orb1419452)

Description	Rabbit polyclonal antibody to GALNTL4.
Species/Host	Rabbit
Reactivity	Bovine, Canine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q6P9A2
Expiration Date	12 months from date of receipt.