

Product Datasheet

CUG-BP1/2 antibody (orb1420010)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to CUG-BP1/2. |
| Species/Host | Rabbit |
| Reactivity | Bovine, Porcine, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC-Fr, IHC-P, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | O95319 |
| Expiration Date | 12 months from date of receipt. |