

Product Datasheet

Guanylyl Cyclase alpha 2 antibody (orb1420440)

Description	Rabbit polyclonal antibody to Guanylyl Cyclase alpha 2.
Species/Host	Rabbit
Reactivity	Bovine, Canine, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P33402
Expiration Date	12 months from date of receipt.