

Product Datasheet

ZNF211 antibody (orb1420851)

Description	Rabbit polyclonal antibody to ZNF211.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13398
Expiration Date	12 months from date of receipt.