

Product Datasheet

MC4 Receptor antibody (orb1421140)

Description	Rabbit polyclonal antibody to MC4 Receptor.
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P32245
Expiration Date	12 months from date of receipt.