

Product Datasheet

MV hemagglutinin glycoprotein antibody (orb1421417)

Description	Rabbit polyclonal antibody to MV hemagglutinin glycoprotein.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08362
Expiration Date	12 months from date of receipt.