

Product Datasheet

Sodium hydrogen exchanger 1 antibody (orb1421432)

Description	Rabbit polyclonal antibody to Sodium hydrogen exchanger 1.
Species/Host	Rabbit
Reactivity	Bovine, Canine, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P19634
Expiration Date	12 months from date of receipt.