

Product Datasheet

Caspase-13 + 4 antibody (orb1421454)

Description Rabbit polyclonal antibody to Caspase-13 + 4.

Species/Host Rabbit

Reactivity Bovine

Conjugation Unconjugated

Tested Applications IF, IHC-Fr, IHC-P, WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Clonality Polyclonal

Uniprot ID [O75601](#)

Expiration Date 12 months from date of receipt.