

Product Datasheet

COVID-19 Nucleocapsid Protein antibody (orb1421466)

Description	Rabbit polyclonal antibody to COVID-19 Nucleocapsid Protein.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P0DTC9
Expiration Date	12 months from date of receipt.