

Product Datasheet

GAPDH antibody (HRP) (orb1421510)

Description	Mouse monoclonal antibody against GAPDH conjugated to HRP.
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Virus
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P04406
Dilution Range	WB: WB (1/2000 - 1/5000)
Expiration Date	12 months from date of receipt.