

Product Datasheet

Beta-tubulin antibody (HRP) (orb1421518)

Description	Mouse monoclonal antibody against Beta-tubulin conjugated to HRP.
Species/Host	Mouse
Reactivity	Human, Insect, Mouse, Rabbit, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IH, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P07437
Dilution Range	WB: WB (1/20000 - 1/600000), IH (1/1000 - 1/2000), IH: WB (1/20000 - 1/600000), IH (1/1000 - 1/2000)
Expiration Date	12 months from date of receipt.