

Product Datasheet

VSV-G-tag antibody (orb1421528)

Description	Mouse monoclonal antibody to VSV-G-tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200), IF: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200), IP: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.