

## Product Datasheet

### V5-tag antibody (orb1421529)

<b>Description</b>	Mouse monoclonal antibody to V5-tag.
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IP, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200), IF: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200), IP: WB (1/2000 - 1/5000), IF/IC (1/200 - 1/500), IP (1/100 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.