

Product Datasheet

TAP-tag antibody (orb1421530)

Description	Mouse monoclonal antibody to TAP-tag.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: WB (1/5000 - 1/10000), IP (1/100 - 1/200), IP: WB (1/5000 - 1/10000), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.