

Product Datasheet

Strep-tag II antibody (orb1421532)

Description	Mouse monoclonal antibody to Strep-tag II.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: WB (1/2000 - 1/5000)
Expiration Date	12 months from date of receipt.