

Product Datasheet

HA-tag antibody (HRP) (orb1421540)

Description	Rabbit polyclonal antibody against HA-tag conjugated to HRP.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: WB (1/1000 - 1/3000)
Expiration Date	12 months from date of receipt.