

Product Datasheet

FLAG-tag antibody (HRP) (orb1421541)

Description	Rabbit polyclonal antibody against FLAG-tag conjugated to HRP.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: WB (1/5000 - 1/10000)
Expiration Date	12 months from date of receipt.