

Product Datasheet

Histone H2B (AcK20) antibody (orb1421678)

Description	Rabbit polyclonal antibody to Histone H2B (AcK20).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96A08, P33778, P62807
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/50 - 1/200), IH: WB (1/500 - 1/1000), IH (1/50 - 1/200), E: WB (1/500 - 1/1000), IH (1/50 - 1/200)
Expiration Date	12 months from date of receipt.