

Product Datasheet

Topoisomerase 2 alpha (pS1106) antibody (orb1421918)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to Topoisomerase 2 alpha (pS1106). |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Porcine, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IH, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P11388 |
| Dilution Range | WB: WB (1/500 - 1/1000), IH (1/50 - 1/200), IH: WB (1/500 - 1/1000), IH (1/50 - 1/200) |
| Expiration Date | 12 months from date of receipt. |