

Product Datasheet

CD127 (pY449) antibody (orb1421921)

Description	Rabbit polyclonal antibody to CD127 (pY449).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P16871
Dilution Range	WB: WB (1/500 - 1/1000), IF/IC (1/100 - 1/500), IF: WB (1/500 - 1/1000), IF/IC (1/100 - 1/500), E: WB (1/500 - 1/1000), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.