

Product Datasheet

Glucocorticoid Receptor (pS226) antibody (orb1422027)

Description	Rabbit polyclonal antibody to Glucocorticoid Receptor (pS226).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IH, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04150
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/50 - 1/100), IH: WB (1/500 - 1/1000), IH (1/50 - 1/100)
Expiration Date	12 months from date of receipt.