

Product Datasheet

MT-ND5 antibody (orb1422382)

Description	Rabbit polyclonal antibody to MT-ND5.
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Virus, Zebrafish
Conjugation	Unconjugated
Tested Applications	IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P03915
Dilution Range	WB: WB (1/500 - 1/1000), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.