

## Product Datasheet

### p38 (pY323) antibody (orb1422611)

<b>Description</b>	Rabbit polyclonal antibody to p38 (pY323).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat, Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IH, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q16539</a>
<b>Dilution Range</b>	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), E: WB (1/500 - 1/1000), IH (1/100 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.