

Product Datasheet

SMAD2 (pS467) antibody (orb1422694)

Description	Rabbit polyclonal antibody to SMAD2 (pS467).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Virus, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IH, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15796
Dilution Range	WB: E (1/5000 - 1/10000), WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IH: E (1/5000 - 1/10000), WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IF: E (1/5000 - 1/10000), WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), E: E (1/5000 - 1/10000), WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.