

## Product Datasheet

### Nogo antibody (orb1422698)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to Nogo.  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rabbit, Rat, Virus   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IH, WB   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Uniprot ID</b>          | <a href="#">Q9NQC3</a>   |
| <b>Dilution Range</b>      | WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200)   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |