

Product Datasheet

HOXC6 antibody (orb1422771)

Description	Rabbit polyclonal antibody to HOXC6.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P09630
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IF: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.